

Wood recliner – SB427A

Seating options designed for today's care environments



Features and benefits

- An extensive selection of fabric and colour configurations, ideal for a range of facilities
- Long-lasting, durable construction, using materials that are both easy to clean and maintain
- Aptly sized to suit a variety of needs across most healthcare settings

Technical Specifications model SB427A

Height	980 - 1,170 mm
Seat height	490 mm
Width	625 mm
Depth	~830 - 1,640 mm
Weight	28 kg
Maximum load*	140 kg
Material	
fabric	coated polycotton blend
foam	high-density fireproof polyurethane
frame	beech plywood
Foam depth	
seat	90 - 135 mm
footrest	35 mm
back - front	10 - 40 mm
back - back	10 mm
Fabrics/colours	
see fabrics and colours sheet	

*as per manufacturer's recommendation

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health**. Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers**.

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN150346-01, 29 February 2016