

Sabina™ 200 sit-to-stand lift

An economical sit-to-stand lift



The Sabina 200 sit-to-stand lift is a basic model sit-to-stand lift, used to assist patients who have difficulty in raising themselves from a seated to a standing position. This can be, for example, from a wheelchair, chair, bed or when using a toilet.

Features and benefits

An extensive assortment of accessories is available

Adjustable height settings accommodate a variety of different patient sizes

Powerful lifting motor ensures the smooth raising and lowering of patients

Base-width adjustment offers greater flexibility of use when lifting from a bed or chair, enhancing patient mobility

The lower leg support and footrest **can easily be removed** for cleaning

Sabina 200 sit-to-stand lift



Powered width expansion



Control box



External charger

Technical Specifications

Maximum Load	160 kg
Weight	42 kg
Protection class	IPX4
Slingbar dimensions	350 mm wide
Base	
height	105 mm
length	1160 mm
width	680 - 1120 mm
expansion	electric
Castors	
front	75 mm double band
rear	75 mm double band
Material	powder coated steel
Emergency lowering	electric
Lifting speed (without load)	20 mm / sec
Lifting interval	795 mm

Sound level	42 dbA
Battery	lead-acid gel
Charger	external
Hand control	wired
Product number	2020020

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

Sabina 200 mobile lift is intended to be used for the lifting and transfer of patients in a variety of care settings.

Class I

Manufacturer: Liko AB, Nedre vägen 100, SE-975 92 Luleå, Sweden

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: 5EN183315-01, 10 February 2017

Enhancing outcomes for patients and their caregivers:

Hill-Rom